

## Vogue Lightweight 2 - Seat Walker



This modern seat walker has a clean and crisp frame design, making it a standout when going for a walk or browsing the shops!

**LIGHTWEIGHT ALUMINIUM BODY:** Improves manoeuvrability and portability.

- Advanced locking breaks
- Unbreakable TPE super soft wheels
- New wider curved backrest
- Minimised storage space with folded width 220 mm
- Locks when folded
- Removable utility bag with 5 kg limit
- Cane holder attachment included
- \*Unit weight includes accessories, bag and backstrap
- \*\*Frame weight is minimal functional weight

✔ Contracted DVA product

✔ Contract approved

### Compliance

ISO 11199-2:2005

## Warranty Information

- 1 Year

| Code        | Display Name                    | Total Length | Total Width | Total Height | Safe Working Load | Unit Weight | Colour    | Seat Width | Seat To Floor Height | Frame Weight |
|-------------|---------------------------------|--------------|-------------|--------------|-------------------|-------------|-----------|------------|----------------------|--------------|
| WAF705310RE | Vogue Lightweight 2 - Red       | 665 mm       | 620 mm      | 785 - 920 mm | 150 kg            | 6.2 kg      | Red       | 465 mm     | 540 mm               | 4.9 kg       |
| WAF705310CH | Vogue Lightweight 2 - Champagne | 665 mm       | 620 mm      | 785 - 920 mm | 150 kg            | 6.2 kg      | Champagne | 465 mm     | 540 mm               | 4.9 kg       |

| Code        | Display Name                  | Total Length | Total Width | Total Height | Safe Working Load | Unit Weight | Colour  | Seat Width | Seat To Floor Height | Frame Weight |
|-------------|-------------------------------|--------------|-------------|--------------|-------------------|-------------|---------|------------|----------------------|--------------|
| WAF705310LE | Vogue Lightweight 2 - Leopard | 665 mm       | 620 mm      | 785 - 920 mm | 150 kg            | 6.2 kg      | Leopard | 465 mm     | 540 mm               | 4.9 kg       |