




RETAIL RANGE



ASPIRE'S RETAIL RANGE

Aspire's new retail box range combines industry leading product with eye catching packaging. The Aspire products come packaged in convenient packets or cardboard boxes with carry handles, they are easy to assemble and come with a minimum one-year warranty.



**Transfer Bench
Plastic Seat & Bench**
Plastic Seat BTS108200
Padded Seat BTS109000
Heavy Duty BTS105300



Bath & Shower Mats
340 x 580mm BTS110200
530 x 530mm BTS109600



Bath Sponge
Plastic Handle BTS111100



Bath Board
BTS113700



Homecare Over Toilet Aid
BTT145025



Homecare Shower Chair
BTS118025



Homecare Shower Stool
BTS119035



Toilet Seat Raiser with Lid
50mm BTT149900
100mm BTT150100
150mm BTT150400



Toilet Seat Raiser without Lid
50mm BTT150200
100mm BTT150300
150mm BTT150400



Male Urinal with Lid
BTT148450



Multi Purpose Cushion
GSS340120



Foam Back Support
GSS340140



Aspire Foam Posture Wedge
GSS340215



Lumbar Roll Cushion
GSS340240



Shoe Horn
DLD269160



Stocking / Sock Aid
DLD269275



Pick Up Reacher
Standard DLR279150
Long DLR279170
Extra Long DLR279230



Heel Protectors
PTA503360



Crutch Accessories
Underarm Cushion
WAA686240
Padded Hand Grip
WAA686250



**Wing Nut
Underarm Crutches**
Youth/Small WAC692503
Adult/Medium WAC692502
Adult Tall/Large WAC692501



**Superlite Adjustable
Seat Walker**
WAF709650CH



Slide Sheet
Small TSE673070
Large TSE673080



Digital Pedal Exerciser
FPP292810



AssistBar
DLG271660



Swivel Cushion
TSE673056



**Aluminium Side Folding
Walking Frame**
WAF710100
Frame, Wheels & Skis - 540mm
WAF710100A



Walking Frame Wheels
Pair WAA686560



Walking Frame Ski's
WAA686550

SUGGESTED FLOOR MERCHANDISED RANGE



Vogue Indoor Walker
 Champagne/Black
 WAF705100CH
 Red/Black
 WAF705100RE



Vogue Lightweight 2 Seat Walker
 Champagne WAF705310CH
 Red WAF705310RE



Vogue Carbon Fibre Seat Walker
 Standard WAF705350
 Tall WAF705360



Vogue Adventure Walker
 WAF705450



Socialite Folding Wheelchair
 Attendant Propelled
 Silver MWS449800
 Red MWS449810



Socialite Folding Wheelchair
 Self Propelled
 Silver MWS449820
 Red MWS449830



Dash Folding Wheelchair
 Attendant Propelled
 Orange MWS449850
 Green MWS449855



Dash Folding Wheelchair
 Self Propelled
 Orange MWS449860
 Green MWS449865



Freedom Excel Superlite Transporter Wheelchair
 Attendant Propelled
 460mm MWS452110
 510mm MWS452120



Freedom Excel Superlite Wheelchair
 Self Propelled
 460mm MWS452170
 510mm MWS452180

