

# Liko™ Sling Sizing Guide

Liko™ Sling	Size	Recommended Weight Range	
Original High Back/Soft Original High Back Mod 20/21, 26 (Soft)	X Small	15 - 20kg	●
	Small	20 - 30kg	●
	Medium Slim	30 - 50kg	●
	Medium	40 - 80kg	●
	Large	70 - 120kg	●
	X Large	110 - 250kg	●
	XXL	200 - 499kg max	○
ComfortSling with Head Support Mod 300/310, 350/360 (High)	Medium	40 - 80kg	●
	Large	70 - 120kg	●
	X Large	100 - 299kg	●
SilhouetteSling with Head Support Mod 22	Medium	45 - 66kg	●
	Large	65 - 100kg	●
Hygiene Sling Mod 40, 41 (Belt), 45/46	X Small	< 30kg	●
	Small	25 - 50kg	●
	Medium	40 - 90kg	●
	Large	80 - 120kg	●
	X Large	100 - 300kg	●
Amputee/High Back Amputee Mod 70, 75 (High Back)	Small	< 40kg	●
	Medium	40 - 60kg	●
	Large	60 - 200kg	●
	X Large (HB Only)	70 - 200kg	●
Ultra Sling™ Mod 6	Short/Regular/Long	Up to 499kg	
Solo High Back Sling (Disposable) Mod 25	Medium Slim	30 - 50kg	●
	Medium	40 - 80kg	●
	Large	70 - 120kg	●
	X Large	110 - 250kg	●
MasterVest™ Mod 64/66	XX Small	< 20kg	●
	X Small	20 - 40kg	●
	Small	40 - 60kg	●
	Medium	60 - 80kg	●
	Large	80 - 100kg	●
	X Large	100 - 130kg	●
	XXL	120 - 300kg max	○

Hygiene Vest  
Mod 50/55 (Teddy)

XX Small	< 20kg	●
X Small	20 - 40kg	●
Small	30 - 60kg	●
Medium	50 - 80kg	●
Large	70 - 100kg	●
X Large	90 - 130kg (300kg max)	●

Lift Pants™ and Ultra Lift Pants™  
Mod 92

X Small	< 30kg	●
Small	20 - 50kg	●
Medium	50 - 85kg	●
Large	70 - 120kg	●
X Large/XXL/XXXL	> 100kg (max 499kg)	○

RepoSheet™/Solo RepoSheet™  
RepoSheet™

Regular	Up to 200kg
Ultra	Up to 500kg

Sabina™ SupportVest™  
Mod 91

Medium	Girth - 142cm max	●
Large	Girth - 173cm max	●
X Large	Girth - 193cm max	●

Sabina™ ComfortVest™  
Mod 95

Small	< 59kg	●
Medium	50 - 79kg	●
Large	70 - 109kg	●
X Large	> 100kg	●

Sabina™ SafetyVest™  
Mod 93

Medium	Girth - 112cm max (80 - 200kg)	●
Large	Girth - 137cm max (70 - 300kg)	●

**Colour Guide:**

- X Small = Red
- Medium = Yellow
- XXL = White
- Small = Orange
- Large = Blue
- Medium Slim = Grey
- XL = Black

See appropriate Liko™ Sling Instruction Guide for complete information.

For Sling Selection Assistance call Aidacare on **1300 133 120** or email [sales@aidacare.com.au](mailto:sales@aidacare.com.au)